



Product Datasheet

Aquaporin 4 antibody (orb688380)

Description	Aquaporin 4 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Synthetic Peptide
Target	AQP4
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P55087
Dilution Range	WB:1:500-1:5000, IHC-P:1:50-1:500, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) HeLa, 2) Mou...



Immunofluorescence analysis of HeLa usin...