


Product Datasheet

MKI67 antibody (orb688382)

Description	MKI67 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthetic Peptide
Target	MKI67
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P46013
Dilution Range	IHC-P:1:100-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse Ki...