



Product Datasheet

GFP antibody (orb688437)


Description	GFP antibody
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	Recombinant Protein
Target	GFP
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa using GFP ...



Immunoprecipitation analysis of HeLa Cel...



Immunoprecipitation analysis of using GF...