



Product Datasheet

Tubulin beta 3 antibody (orb688448)


Description	Tubulin beta 3 antibody
Species/Host	Mouse
Reactivity	Canine, Gallus, Hamster, Human, Monkey, Mouse, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Synthetic Peptide
Target	TUBB3
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	Q13509
Dilution Range	WB:1:500-1:5000, IHC-P:1:50-1:500, IF:1:50-1:200, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) HeLa, 2) Mou...



Immunoprecipitation analysis of Mouse Br...



Immunohistochemical staining of Human CO...