

Product Datasheet

TUBB antibody (HRP) (orb688489)

Description	TUBB antibody (HRP)
Species/Host	Mouse
Reactivity	Canine, Gallus, Hamster, Human, Insect, Monkey, Mouse, Rabbit, Rat, Sheep, Yeast
Conjugation	HRP
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Protein
Target	TUBB
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P07437
Dilution Range	WB:1:500-1:5000, IHC-P:1:50-1:500
Expiration Date	12 months from date of receipt.