

Product Datasheet

Aquaporin 4 antibody (orb688514)

Description	Aquaporin 4 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Protein
Target	AQP4
Preservatives	Mouse IgG1 in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P55087
Dilution Range	WB:1:500-1:5000, IHC-P:1:50-1:500, IF:1:50-1:200
Expiration Date	12 months from date of receipt.