



## Product Datasheet

### BRAP antibody (orb688531)

|                            |                                                                                                        |
|----------------------------|--------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | BRAP antibody                                                                                          |
| <b>Species/Host</b>        | Mouse                                                                                                  |
| <b>Reactivity</b>          | Human                                                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                                                           |
| <b>Tested Applications</b> | ELISA, WB                                                                                              |
| <b>Immunogen</b>           | Purified recombinant human BRAP protein fragments expressed in E.coli.                                 |
| <b>Target</b>              | BRAP                                                                                                   |
| <b>Preservatives</b>       | Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol. |
| <b>Form/Appearance</b>     | Liquid                                                                                                 |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                          |
| <b>Note</b>                | For research use only                                                                                  |
| <b>Isotype</b>             | IgG1                                                                                                   |
| <b>Clonality</b>           | Monoclonal                                                                                             |
| <b>Uniprot ID</b>          | <a href="#">Q7Z569</a>                                                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                        |



Western blot analysis of CHO-K1 cells us...



Western blot analysis of BRAP in using B...