

## Product Datasheet

### CDK1 antibody (orb688541)

|                            |                                                                                                                             |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | CDK1 antibody                                                                                                               |
| <b>Species/Host</b>        | Mouse                                                                                                                       |
| <b>Reactivity</b>          | Canine, Human, Mouse, Rat                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                                                |
| <b>Tested Applications</b> | ELISA, ICC, WB                                                                                                              |
| <b>Immunogen</b>           | Recombinant fragment, corresponding to amino acids 8-298aa of human CDC2.                                                   |
| <b>Target</b>              | CDK1                                                                                                                        |
| <b>Preservatives</b>       | Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol |
| <b>Form/Appearance</b>     | liquid                                                                                                                      |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                               |
| <b>Note</b>                | For research use only                                                                                                       |
| <b>Isotype</b>             | IgG2b                                                                                                                       |
| <b>Clonality</b>           | Monoclonal                                                                                                                  |
| <b>Uniprot ID</b>          | <b>P06493</b>                                                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                             |