

Product Datasheet

FOS antibody (orb688546)

Description	FOS antibody
Species/Host	Mouse
Reactivity	Human, Monkey, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified recombinant human FOS protein fragments expressed in E.coli.
Target	FOS
Preservatives	Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P01100
Expiration Date	12 months from date of receipt.