

Product Datasheet

CHRD1 antibody (orb688550)

Description	CHRD1 antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified recombinant human CHRD1 protein fragments expressed in E.coli.
Target	CHRD1
Preservatives	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q9BU40
Expiration Date	12 months from date of receipt.