

Product Datasheet

CRK antibody (orb688554)

Description	CRK antibody
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified recombinant human CRK protein fragments expressed in E.coli.
Target	CRK
Preservatives	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P46108
Expiration Date	12 months from date of receipt.