

Product Datasheet

DLAT antibody (orb688560)

Description	DLAT antibody
Species/Host	Mouse
Reactivity	Human, Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IP, WB
Immunogen	Purified recombinant human DLAT protein fragments expressed in E.coli.
Target	DLAT
Preservatives	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P10515
Dilution Range	WB:1:500-1:5000, ICC:1:50-1:200, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.