

Product Datasheet

EGFP/EYFP antibody (orb688565)

Description	EGFP/EYFP antibody
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified recombinant EYFP.tag full length expressed in E.coli.
Target	EGFP/EYFP
Preservatives	Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	C8CHS1
Expiration Date	12 months from date of receipt.