

Product Datasheet

FAS antibody (orb688567)

Description	FAS antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant fragment, corresponding to amino acids 211aa-329aa of human CD95(C-terminus).
Target	FAS
Preservatives	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P25445
Expiration Date	12 months from date of receipt.