

Product Datasheet

FYB antibody (orb688571)

Description	FYB antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	Purified recombinant human FYB protein fragments expressed in E.coli.
Target	FYB
Preservatives	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	O15117
Expiration Date	12 months from date of receipt.