

Product Datasheet

GAPDH antibody (orb688572)

Description	GAPDH antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant fragment, corresponding to amino acids 72aa-248aa of human GAPDH (Internal).
Target	GAPDH
Preservatives	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P04406
Expiration Date	12 months from date of receipt.