

Product Datasheet

IL18BP antibody (orb688588)

Description	IL18BP antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified recombinant human IL18BP protein fragments expressed in E.coli.
Target	IL18BP
Preservatives	Ascites
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	O95998
Expiration Date	12 months from date of receipt.