

Product Datasheet

CCAR2 antibody (orb688596)

Description	CCAR2 antibody
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant fragment, corresponding to amino acids 1aa-265aa of human KIAA1967(N-terminus).
Target	CCAR2
Preservatives	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q8N163
Expiration Date	12 months from date of receipt.