

Product Datasheet

MAP2 antibody (orb688599)

Description	MAP2 antibody
Species/Host	Mouse
Reactivity	Bovine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	Recombinant fragment, corresponding to amino acids 1636aa-1855aa of human MAP2(C-terminus).
Target	MAP2
Preservatives	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P11137
Expiration Date	12 months from date of receipt.