

Product Datasheet

CD4 antibody (orb688638)

Description	CD4 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Recombinant Protein
Target	CD4
Preservatives	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P01730
Expiration Date	12 months from date of receipt.