

Product Datasheet

CD4/CD8 antibody (FITC/PE) (orb688640)

Description	CD4/CD8 antibody (FITC/PE)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Immunogen	Recombinant Protein
Target	CD4/CD8
Preservatives	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) rAlbumin
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2b/IgG1
Clonality	Monoclonal
Uniprot ID	P01730
Expiration Date	12 months from date of receipt.