

Product Datasheet

IgG1 antibody (FITC) (orb688663)

Description	IgG1 antibody (FITC)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Recombinant Protein
Target	IgG1
Preservatives	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) rAlbumin
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.