

Product Datasheet

PPP5C antibody (orb688684)

Description	PPP5C antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified recombinant human PPP5C protein fragments expressed in E.coli.
Target	PPP5C
Preservatives	Ascites
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P53041
Expiration Date	12 months from date of receipt.