

Product Datasheet

CD133 antibody (orb688686)

Description	CD133 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified recombinant human PROM1 protein fragments expressed in E.coli.
Target	PROM1
Preservatives	Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide,0.1mg/mlrAlbumin and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	O43490
Expiration Date	12 months from date of receipt.