

## Product Datasheet

### RANBP9 antibody (orb688690)

|                            |                                                                                                                             |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | RANBP9 antibody                                                                                                             |
| <b>Species/Host</b>        | Mouse                                                                                                                       |
| <b>Reactivity</b>          | Human                                                                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                                                                |
| <b>Tested Applications</b> | ELISA, ICC, IP, WB                                                                                                          |
| <b>Immunogen</b>           | Purified recombinant human RANBP9 protein fragments expressed in E.coli.                                                    |
| <b>Target</b>              | RANBP9                                                                                                                      |
| <b>Preservatives</b>       | Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol |
| <b>Form/Appearance</b>     | liquid                                                                                                                      |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                               |
| <b>Note</b>                | For research use only                                                                                                       |
| <b>Isotype</b>             | IgG1                                                                                                                        |
| <b>Clonality</b>           | Monoclonal                                                                                                                  |
| <b>Uniprot ID</b>          | <a href="#">Q96S59</a>                                                                                                      |
| <b>Dilution Range</b>      | WB:1:500-1:5000, ICC:1:50-1:200, IP:1:200-1:2000                                                                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                             |