

Product Datasheet

CD8A antibody (FITC) (orb688693)

Description	CD8A antibody (FITC)
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Recombinant Protein
Target	CD8A
Preservatives	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) rAlbumin
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P01732
Expiration Date	12 months from date of receipt.