

Product Datasheet

TBL1Y antibody (orb688707)

Description	TBL1Y antibody
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	Recombinant full length protein (522aa) of human TBL1Y.
Target	TBL1Y
Preservatives	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	Q9BQ87
Expiration Date	12 months from date of receipt.