

Product Datasheet

GTF2H1 antibody (orb688711)

Description	GTF2H1 antibody
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	Recombinant fragment, corresponding to amino acids 1aa-279aa of human GTF2H1 (N-terminus).
Target	GTF2H1
Preservatives	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P32780
Expiration Date	12 months from date of receipt.