



Product Datasheet

ARNT antibody (orb688728)

Description	ARNT antibody
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Protein
Target	ARNT
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P27540
Dilution Range	WB:1:500-1:5000, IHC-P:1:50-1:500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse Brain tissue.



Western blot analysis of Mouse Brain tissue.