



Product Datasheet

PPARD antibody (orb688731)

Description	PPARD antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Protein of PPAR Delta of PPAR Delta
Target	PPARD
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	Q03181
Dilution Range	IHC-P:1:100-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse Brain tissue



Immunohistochemical staining of Human Brain tissue