


Product Datasheet

TUBE1 antibody (orb688745)

Description	TUBE1 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic Peptide of Epsilon Tubulin
Target	TUBE1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	Q9UJT0
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) MCF7Cell usi...