



Product Datasheet

CBX3 antibody (orb688758)


Description	CBX3 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Protein of HP-1γ
Target	CBX3
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	Q13185
Dilution Range	WB:1:500-1:5000, IHC-P:1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) HeLa Cell Ly...



Immunohistochemical staining of Human Co...



Immunohistochemical staining of human pl...