



Product Datasheet

KRT14 antibody (orb688848)

Description	KRT14 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthesized peptide derived from human Cytokeratin Pan
Target	KRT14/KRT16/KRT5/KRT6A/KRT8
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2b, Kappa
Clonality	Monoclonal
Uniprot ID	P02533 , P02538 , P08779 , P05787 , P13647
Dilution Range	IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human sk...



Immunohistochemical staining of Rat Test...