



Product Datasheet

CD31 antibody (orb688857)


Description	CD31 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Immunogen	Recombinant Human CD31 protein (28-315AA)
Target	CD31
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P16284
Dilution Range	WB:1:500-1:5000, IHC-P:1:50-1:500, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of HL-60 cells ...



Western blot analysis of Human Brain Tis...



Immunohistochemical staining of human to...