



Product Datasheet

CDH1 antibody (orb688858)


Description	CDH1 antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Immunogen	Recombinant Human Cadherin-1 protein (155-707AA)
Target	CDH1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P12830
Dilution Range	WB:1:500-1:5000, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of MCF-7 cells ...



Western blot analysis of Hela Cell Lysat...



Immunofluorescence analysis of Hela cell...