



Product Datasheet

OCT4 antibody (orb688861)


Description	OCT4 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Immunogen	Recombinant Human POU domain, class 5, transcription factor 1 protein (1-360AA)
Target	OCT4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	Q01860
Dilution Range	WB:1:500-1:10000, IHC-P:1:50-1:200, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse brain tis...



Immunohistochemical staining of human lu...



Immunohistochemical staining of human te...