



Product Datasheet

CD14 antibody (orb688867)


Description	CD14 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Immunogen	Recombinant Human Monocyte differentiation antigen CD14 protein (20-345AA)
Target	CD14
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P08571
Dilution Range	WB:1:1000-1:5000, IHC-P:1:50-1:200, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
RAW264.7
cel...



Western blot
analysis of
Human Colon
Car...



Western blot
analysis of
Hela whole
cell...