



Product Datasheet

HSPA8 antibody (orb688870)


Description	HSPA8 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, IP, WB
Immunogen	Recombinant Human Heat shock cognate 71 kDa protein (2-646AA)
Target	HSPA8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P11142
Dilution Range	WB:1:10000-1:256000, IHC-P:1:100-1:500, FC:1:100-1:300, IP:1 µl-4 µl
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pr...



Immunohistochemical staining of human pr...



Immunoprecipitation analysis of HeLa who...