



Product Datasheet

HSPA8 antibody (orb688871)


Description	HSPA8 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Immunogen	Recombinant Human Heat shock cognate 71 kDa protein (2-646AA)
Target	HSPA8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P11142
Dilution Range	WB:1:5000-1:32000, IHC-P:1:100-1:500, IF:1:100-1:300, FC:1:100-1:300
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of MCF-7 cells ...



Western blot analysis of Human Colon Car...



Western blot analysis of HeLa whole cell...