



Product Datasheet

ENO1 antibody (orb688875)


Description	ENO1 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, IP, WB
Immunogen	Recombinant Human Alpha-enolase protein (2-434AA)
Target	ENO1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P06733
Dilution Range	WB: 1:5000-1:320000, IHC-P:1:200-1:500, IF:1:50-1:200, FC:1:100-1:300, IP:1 µl-4 µl
Expiration Date	12 months from date of receipt.



Immunoprecipitation analysis of HepG2 whole cel...



Western blot analysis of Human Colon Carcinoma cell lines



Western blot analysis of MCF-7 whole cell lysate