



Product Datasheet

CD45 antibody (orb688878)


Description	CD45 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Immunogen	Recombinant Human Receptor-type tyrosine-protein phosphatase C protein (24-575AA)
Target	CD45
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P08575
Dilution Range	WB: 1:2000-1:32000, IHC-P: 1:200-1:600, IF: 1:100-1:300, FC: 1:200-1:600
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Raji cells u...



Western blot
analysis of
Human Colon
Car...



Western blot
analysis of
using CD45
anti...