

Product Datasheet

GAPDH antibody (orb688921)

Description	GAPDH antibody	Image not found or type unknown
Species/Host	Mouse	Western blot
Reactivity	Human, Mouse, Rabbit, Rat	analysis of
Conjugation	Unconjugated	Human Colon
Tested Applications	ELISA, IF, IHC, IP, WB	Car...
Immunogen	Recombinant Human GAPDH protein (2-335AA)	
Target	GAPDH	
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol	
Form/Appearance	Liquid	
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.	
Note	For research use only	
Isotype	IgG1	
Clonality	Monoclonal	
Dilution Range	WB: 1:5000-1:80000, IHC-P: 1:200-1:500, IF: 1:50-1:100, IP: 2µl-8µl	
Expiration Date	12 months from date of receipt.	