

Product Datasheet

Goat Anti-Rabbit IgG(H+L) antibody (Biotin) (orb688927)

Description	Goat Anti-Rabbit IgG(H+L) antibody (Biotin)
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Biotin
Tested Applications	ELISA, FC, IF, IHC, IP, WB
Immunogen	Rabbit IgG(H+L)
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:5000-20000, IHC-P:1:200-2000, IF:1:200-2000, IP:1:200-1:2000, FCM:1:200-2000, ELISA: 1:5000-20000
Expiration Date	12 months from date of receipt.