



Product Datasheet

TRAF2 antibody (orb688997)


Description	Rabbit monoclonal antibody to TRAF2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, IP, WB
Immunogen	A synthesized peptide derived from human TRAF2
Target	TRAF2
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	Q12933
Dilution Range	WB:1:500-1:5000, IHC-P:1:50-1:200, FC:1:20-1:200, IP:1:200-1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human kidney tissue.



Immunohistochemical staining of human cells.



Flow cytometric analysis of HeLa Cells.