

Product Datasheet

Neutrophil Elastase antibody (orb689028)

Description	Rabbit monoclonal antibody to Neutrophil Elastase
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC
Immunogen	A synthesized peptide derived from human Neutrophil Elastase
Target	ELANE
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P08246
Dilution Range	IHC-P:1:50-1:200, FC:1:20-1:200
Expiration Date	12 months from date of receipt.

