



## Product Datasheet

### KLF4 antibody (orb689029)


|                            |                                                                                                   |
|----------------------------|---------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit monoclonal antibody to KLF4                                                                |
| <b>Species/Host</b>        | Rabbit                                                                                            |
| <b>Reactivity</b>          | Human                                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                                      |
| <b>Tested Applications</b> | ELISA, IF, IHC, IP, WB                                                                            |
| <b>Immunogen</b>           | A synthesized peptide derived from human KLF4                                                     |
| <b>Target</b>              | KLF4                                                                                              |
| <b>Preservatives</b>       | Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| <b>Form/Appearance</b>     | Liquid                                                                                            |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                     |
| <b>Note</b>                | For research use only                                                                             |
| <b>Isotype</b>             | Rabbit IgG                                                                                        |
| <b>Clonality</b>           | Monoclonal                                                                                        |
| <b>Uniprot ID</b>          | <b>O43474</b>                                                                                     |
| <b>Dilution Range</b>      | WB:1:500-1:5000, IHC-P:1:50-1:200, IF:1:20-1:200, IP:1:200-1:1000                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                   |



Immunoprecipitation analysis of HeLa whole cell lysate...



Western blot analysis of HeLa whole cell lysate...



Immunohistochemical staining of human skin tissue...