



Product Datasheet

VCP antibody (orb689052)


Description	Rabbit monoclonal antibody to VCP
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	A synthesized peptide derived from human VCP
Target	VCP
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P55072
Dilution Range	WB:1:500-1:5000, IHC-P:1:50-1:200, IF:1:20-1:200, IP:1:200-1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human glioma tissue.



Immunofluorescence analysis of SY5Y cells.



Immunoprecipitation analysis of U251 cells.