



Product Datasheet

PABPN1 antibody (orb689055)


Description	Rabbit monoclonal antibody to PABPN1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Immunogen	A synthesized peptide derived from human PABPN1
Target	PABPN1
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	Q86U42
Dilution Range	WB:1:500-1:5000, IHC-P:1:50-1:200, IF:1:20-1:200, FC:1:20-1:200
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of hela cells
u...



Western blot analysis of 293 whole cell ...



Immunohistochemical staining of human te...