



Product Datasheet

SDHB antibody (orb689064)


Description	Rabbit monoclonal antibody to SDHB
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, WB
Immunogen	A synthesized peptide derived from human SDHB
Target	SDHB
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P21912
Dilution Range	WB:1:500-1:5000, IHC-P:1:50-1:200, FC:1:20-1:200
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of 293 cells us...



Western blot analysis of 293 whole cell ...



Immunohistochemical staining of human he...