



Product Datasheet

BUB1B antibody (orb689068)

Description	Rabbit monoclonal antibody to BUB1B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Immunogen	A synthesized peptide derived from human BubR1
Target	BUB1B
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	O60566
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 whole cell...



Flow cytometric analysis of HeLa cells u...